## CT Whole Body Phantom **PBU-60**

Product No. PH-2B

Frog leg position

### Full internal organs for CT scanning with appropriate HU numbers.


# A unique, life size whole body phantom for CT provides a variety of educational application as well as visual evaluation in finding out optimal scanning conditions. The phantom can also be used for plain X-ray, showing life-like images. No metal parts or liquid structure are used.



Whole Body CT







Main joints have close-to human articulation, allowing various positioning for training.




rotate through a full 360 degrees in the sagittal plane, approx. 180 degrees to side-ways.





rotate forward up to approx. 90 degrees, then abduct up to 45 degrees each.







**Incorporated Anatomy** 

### **Internal Organs** Head

Synthetic skull

Cervical vertebrae Mesencephalon Cerebellum Cerebral ventricles Eye balls Arteries with contrast medium (left half only)
























1.089



Trunk

Clavicles Coxal bones Femurs Lungs with pulmonary vessels (up to third bifurcations) Trachea (up to fourth bifurcations) Liver with portal and hepatic veins Pancereas Kidnevs Gallbladder Urinary bladder Sigmoid Colon



Materials and features:

#### Soft tissue (SZ-50) material

Hounsfield number (Approximation)

Soft tissue and organs: Urethane base resin (SZ-50) Synthetic bones: Joint attachments: Epoxy, urethane with carbon fiber Screws: Poly carbonate

Material	Density g/cm3	Effective Atomic Number	Electron Density x10~23e/g	Elemental Composition (wt%)			
				Н	С	N	0
Water	1.000	7.417	3.343	11.19			88.81
SZ-50	1.061	6.14	3.258	8.41	72.25	4.61	14.73
		Soft tissue		Liver		Kidneys	

1.061

KYOTO KAGAKU

MedicalQC

深圳为尔康科技有限公司 联系人:曾祥满 手机:13632925349

QQ: 274798107 电话: 0755-28896837 地址: 深圳市龙岗区沙平北路111号608A

网址:www.medicalQC.com 邮箱:szchina1718@163.com

1.075